

## GOI-1 Picosecond Gated Optical Image Intensifier

### Single Channel GOI



**Applications** Plasma Diagnostics, time dependent interferometry, etc.

### Summary

- Solid sate drive electronics
- 18 mm diameter input and output format
- Minimum exposure time <100ps typically ~80ps</li>
- Four operational modes
- Computer control
- Trigger delay adjustment over 50ns range
- Cathodes available through visible and near infra red.
- Electronics in 19 inch rack mount up to 3 metres from optical head
- · Custom configurations available

The GOI-1 is a turnkey system comprised of solid-state electronics, an18mm image intensifier ready to be fitted with the users readout equipment and input optics. This latest model feature high spatial resolution at the fastest gate times. High-speed pulsers deliver typical minimum exposures of 80ps FWHM. The design of the electronics allows performance with very low jitter and a small trigger delay. The design of the system allows easy alteration of gate times and easy upgrades to accommodate additional GOI heads with very low channel to channel jitter. With 80ps gate widths, users may extract the most relevant information from their experiments and also have better spatial resolution. Users can precisely control the trigger delay from 0 to 50ns via a computer controlled passive delay line which introduces no further jitter.

The camera may be operated in four different modes: DC on, slow gate, medium gate, and fast gate. The DC mode allows easy focusing and set-up.

The standard offerings for the input window include: glass, fiber optic, or quartz for an extended blue response. Users can select from S20, S25, or S1 photocathodes, providing a solution to demands for imaging in the visible and NIR.

See www.kentech.co.uk



### GOI-1 Picosecond Gated Optical Image Intensifier

SPECIFICATIONS

Number of channels 1 x 18mm wafer micorchannel plate Gen-II tube.
Tube sensitivity Tubes can be selected with S1, S20 S25 type response.

Input window Quartz as standard with S20, others on glass.

Optical output P43 phosphor on fibre optic window

Gate modes Fast, Medium, Slow and DC [cathode on for 5 seconds from software trigger]

Optical gate widths Fast mode nominally: <100, 100, 120 ps

Medium mode: 250, 500, 1000, 2000, 3000, 4000, 5000 ps

Slow mode: 100ns to 1ms ~54 ns in fast mode.

Jitter <20ps Standard Deviation, typically ~4ps.

Gate time trigger delay range ~50 ns in ~25ps steps.

Gain control MCP voltage (typically 260 to 900 volts) is mapped to linear gain settings

0 through  $\bar{1}000$ 

Maximum repetition rate 100Hz.

External trigger requirements 5 volts into  $50\Omega$  rising in < 5ns for optimum jitter.

or 500 volts into  $50\Omega$  with <1ns rise [used when synchronising several channels]

**Controls** 

Control RS232 serial port or Ethernet port. [USB to special order].
Commands Simple structured text commands formatted for easy parsing.

**Power** Switches AC power in the pulser.

Input power range 100 to 240 VAC 50 to 60Hz at <100 watts.

**Indicators** 

Power Shows that AC power is applied and the unit is switched on.

Triggered Illuminates when the unit is triggered. Mode 1 of 4 LEDs will be illuminated.

Overload Illuminates if too much current is drawn from the the high voltage power supply

Dimensions (control unit)

width 19 inch rack mount, 482 mm over handles,

depth 500mm into 19 inch rack, 540 mm over handles, plus rear cover adding 45 mm.

height 128.5 mm (3U)

weight 7.5 kg

**Dimensions** (Optical head)

 height
 132mm

 width
 102mm

 depth
 ~35.5mm

umbilical dimameter nominally 20mm, length 3 metres.

Connectors

Power inlet IEC Trigger input BNC

Head Combination of 1 x Lemo, 2 x SMA and 1 x GES HV connectors

**Environmental issues** 

Storage The unit should be stored in an environment within the following parameters:

Temperature between 10° C and 40° C Humidity <60% non-condensing

Pressure 40 to 120 kPa
Gas type Air or nitrogen

Use The unit should be used in an environment within the following parameters:


Temperature between 10° C and 30° C
Humidity < 60% on condensing
Pressure 80 to 120 kPa
Gas type Air or nitrogen

See www.kentech.co.uk





# GOI-1 Picosecond Gated Optical Image Intensifier

Shortest Gate 85.7ps FWHM = 76ps after correction for laser diode pulse length



A medium range gate Gate ~500ps FWHM





See www.kentech.co.uk

